


CT Scan Protocol

Knee

The CT scan quality is critical to the production of accurate personalized implants and patient-specific guides. Deviations from this protocol may result in an unusable scan and delay of surgery. Please contact Meticuly team for further clarification.

Scanning Parameters

Region of interest	Distal 1/3 of femur through proximal 1/2 of tibia with entire patella (approx. 25 cm above and below knee)	
Body side	Both left and right legs	
Matrix size	512 x 512	
Voxel size	0.5 – 1 mm	
Slice thickness	maximum 1 mm	
Feed per rotation	maximum 1 mm	
Pitch	1 or less	
Reconstructed slice increment	maximum 1 mm	
Reconstruction algorithm	Bone	
Export File	DICOM	
File Format	Uncompressed standard	

CT Scanning Instruction

- Helical (spiral) scanning mode is preferred for CT image acquisition. A conventional CT can be used if minimum requirements stated above could be fulfilled.
- Both left and right sides should be scanned with approximately the same setting
- Position the patient to maximize comfort and minimize motion. If possible, position the patient as follows: supine, feet first, patellae pointing forward and the knees in maximal extension, toes pointing straight up.
- Use the smallest field of view possible to capture the whole regions of the required bones. Capturing all soft tissue is unnecessary, only the bony regions are of interest.
- Images scanned with no gantry tilt and no oblique reconstruction (i.e. use only primary axial images). No reformatting into coronal or sagittal planes.
- All slices must have the same field of view, reconstruction center, and table height.
- Scan with the same slice spacing, less than or equal to the slice thickness.

Data Transfer

- Provide the complete data set of raw/original DICOM images to the surgeon
- Do not erase patient name and ID. Data will be anonymized by Meticuly on receipt of the data, after cross-check with prescription of the surgeon to ensure the images of the right patient are provided.